Structural basis of AAV affinity chromatography: from design to cryo-EM validation

Nathaniel E Clark¹, Mario Mietzsch², Robert McKenna², Kelley Kearns¹, Thomas C Scanlon¹

nclark@repligen.com

Introduction

To increase patient access to AAV therapies, it is necessary to decrease manufacturing costs. To this end, we developed a suite of alkali-tolerant affinity resins for AAV serotypes 2, 5, 6, 8, and 9. The resulting affinity resins offer best-in-class performance and process economics as they maintain performance through 20+ clean-in-place cycles.

To validate that design and specificity of our AVIPure[®] AAV-affinity ligands we determined cryo-EM structures of the AAV2, AAV8, and AAV9 ligands in complex with capsids of the appropriate serotype.

The structures map the footprint of AVIPure ligands on the capsids, enabling rational pairing of affinity resin to an engineered capsid. If capsid engineered does not substitute amino acids in AVIPure footprint, successful purification with AVIPure resin is expected.

AVIPure scaffolds can bind virtually any biological epitope

AVIPure ligands are: 1) Stable in NaOH 2) Small proteins and peptides 4) Diverse set of >49 scaffolds 3) Manufacturable at scale

AVIPure[®] Affinity Resin portfolio and pipeline







Day 1	
Week 2	
Week 3	
Week 4	
Week 6	

\checkmark				
		BIND		
	4000 - 3000 - (NPW) 2000 - 2000 - 1000 - 0 -		6	12
			•	





Learn more at **repligen.com**



¹Repligen Corporation, 16 Cavendish Ct Lebanon NH 03766 ²University of Florida, Gainesville FL 32610

REPLIGEN

© 2023 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners