



PATsmart™ REBEL® System

Analyte Panel Performance Guide

26 Nov 2025

PATsmart™ REBEL® System Analysis—Calibration Range

SMA v2: Spent Media Analysis Kit, version 2

- As with all analytical assays, it is vital to ensure that analyte concentrations are within the appropriate dynamic range.
- Analyte concentrations below or above the calibration range will exhibit decreased reproducibility and accuracy.
- REBEL software has a built-in dilution guide— Range Finding Experiment Tool (RFE)—providing automated assessment of optimal dilution factors for the most inclusive analysis.


| SMA V2 Analytes | Abbreviation | Lower Limit of Quantitation (LLOQ) | Upper Limit of Quantitation (ULOQ) |
|---------------------|--------------|------------------------------------|------------------------------------|
| Alanyl-Glutamine | AQ | 5 µM (0.005 mM) | 100 µM (0.100 mM) |
| Alanine | Ala | | |
| Arginine | Arg | | |
| Asparagine | Asn | | |
| Aspartic acid | Asp | | |
| Betaine | Betaine | | |
| Choline | Choline | | |
| Citrulline | Cit | | |
| Glycine | Gly | | |
| Histidine | His | | |
| Hydroxyproline | Hyp | | |
| Isoleucine | Ile | | |
| Leucine | Leu | | |
| Lysine | Lys | | |
| l-Methylhistidine | lMH | | |
| Methionine | Met | | |
| Phenylalanine | Phe | | |
| Proline | Pro | | |
| Serine | Ser | | |
| Threonine | Thr | | |
| Tryptophan | Trp | | |
| Tyrosine | Tyr | | |
| Valine | Val | | |
| Glutamine | Gln | 5 µM (0.005 mM) | 75 µM (0.075 mM) |
| Glutamic acid | Glu | 5 µM (0.005 mM) | 50 µM (0.050 mM) |
| Cystine | Cystine | | |
| β-Alanine | βAla | 5 µM (0.005 mM) | 25 µM (0.025 mM) |
| Thiamine | B1 | | |
| Pyridoxine | B6-OH | | |
| Pyridoxal | B6-Oxo | | |
| Nicotinamide | NAM | | |
| β-Aminobutyric acid | GABA | | |
| Sarcosine | Sarcosine | | |

REBEL Error Reporting

Understanding standard error on the REBEL

- REBEL reports all results in mM concentration accompanied by the associated standard error.
- Standard error and relative standard error (RSE) are statistical propagations of error for the system as a whole, based on the QC and individual sample replicate.
- Not to be confused with Coefficient of Variation (%CV) or Standard Deviation.
- The greater the RSE, the greater the uncertainty in the results .
 - a) Most common cause for large RSE is extrapolation outside of the calibrated range, or at the extremes of the curve.
 - b) Acceptable RSE are established by each user group, but generally <15% is considered good.


Relative standard error (RSE) can be interpreted as the expected variability in estimated concentrations over many replicates—samples are recommended to be run in triplicate.



REBEL Analysis—Calibration Range

Determining whether analyte data is in/out of range

- When preparing samples, choosing the optimal dilution factor prior to analysis is vital to ensuring precise and accurate results. Range Finding Experiment (RFE) tool helps optimize dilution factor for unknown samples.
- The proper dilution should bring samples in range so results fall within the standard calibration curve post-dilution or risk lower confidence/higher relative standard error (RSE) in resulting data.
- All samples require a minimum of at least 10x dilution in REBEL diluent.



Analytes frequently in low levels in cell culture media

REBEL analyte panel considerations

- Vitamins and other amines in the analyte panel may present challenges based on their conventionally lower concentrations in most media, and the smaller dynamic range within which the analytes are calibrated for analysis.
- Further optimization of dilution factor distinct from AAs may be required—i.e., a lower dilution factor is often necessary to accurately quantify these analytes.

| Analyte full name | Abbreviation |
|---------------------|--------------|
| Vitamins | |
| Choline | Choline |
| Pyridoxine | B6-OH |
| Pyridoxal | B6-Oxo |
| Nicotinamide | NAM |
| Thiamine | B1 |
| Amines | |
| γ-aminobutyric acid | GABA |
| 1-Methylhistidine* | 1MH |
| Citrulline | Cit |
| Hydroxyproline | Hyp |
| Sarcosine | Sarcosine |

*1-methylhistidine cannot be distinguished from 3-methylhistidine

Analytes with specific considerations in cell culture media

REBEL analyte panel considerations

- For optimal performance of all analytes, please use standards and reagents within expiration.
 - > Posted Expiration—BGE/Diluent expire 1 year after manufacture; Standards expire 3 months after manufacture.
 - > Use Expiration—BGE/Diluent expire 30 days after opening; Standards expire 14 days after initial use.

| Analyte full name | Abbreviation | Comment |
|-------------------|--------------|--------------------------------------|
| β-Alanine | βAla | Measurements may be less precise |
| Arginine | Arg | Measurements may be less precise |
| Cystine | Cystine | May be prone to degradation |
| Glutamine | Gln | Possible co-migration with Glu |
| Glutamic acid | Glu | Possible co-migration with Gln |
| Glycine | Gly | Measurements may be less precise |
| Serine | Ser | Measurements may be less precise |
| Valine | Val | Possible co-migration with Sarcosine |

- “Measurements may be less precise”—under certain conditions or in some media, the measurements resulting for this analyte may be less reproducible than others.
- “May be prone to degradation”—analyte is known to be more labile under certain conditions.
- “Possible co-migration”—analyte known to encounter detection interferences with other analytes, especially at higher concentrations.

Summary

REBEL performance is expected to be accurate and precise, while there are notable analytes that may prove more challenging based on observations, this information is good to keep in mind when analyzing REBEL data or preparing dilutions.

Routine Amino Acids

| Analyte full name | Abbreviation |
|---------------------------|--------------|
| Alanyl-Glutamine | AQ |
| Alanine | Ala |
| Asparagine | Asn |
| Aspartic acid | Asp |
| Glutamine + Glutamic Acid | Glx |
| Histidine | His |
| Isoleucine | Ile |
| Leucine | Leu |
| Lysine | Lys |
| Methionine | Met |
| Phenylalanine | Phe |
| Proline | Pro |
| Threonine | Thr |
| Tryptophan | Trp |
| Tyrosine | Tyr |

Vitamins + Amines¹

| Analyte full name | Abbreviation |
|--------------------------------|--------------|
| Vitamins | |
| Choline | Choline |
| Pyridoxine | B6-OH |
| Pyridoxal | B6-Oxo |
| Nicotinamide | NAM |
| Thiamine | B1 |
| Amines | |
| γ-aminobutyric acid | GABA |
| 1-Methylhistidine ² | 1MH |
| Citrulline | Cit |
| Hydroxyproline | Hyp |
| Sarcosine | Sarcosine |

¹ These analytes may present challenges based on their lower concentrations in media and smaller dynamic range for analysis.
² 1-methylhistidine cannot be distinguished from 3-methylhistidine.

Analytes with specific considerations

| Analyte full name | Abbreviation | Comment |
|-------------------|--------------|--------------------------------------|
| β-Alanine | βAla | Measurements may be less precise |
| Arginine | Arg | Measurements may be less precise |
| Cystine | Cystine | May be prone to degradation |
| Glutamine | Gln | Possible co-migration with Glu |
| Glutamic acid | Glu | Possible co-migration with Gln |
| Glycine | Gly | Measurements may be less precise |
| Serine | Ser | Measurements may be less precise |
| Valine | Val | Possible co-migration with Sarcosine |

Expirations and LOQ

| Reagent | Expiration |
|---------------|--|
| Standards | 14 days after initial use (or 3 months from manufacture – posted date) |
| BGE / Diluent | 30 days after opening (or 1 year from manufacture – posted date) |

| LLOQ/ULOQ (mM) | Analytes |
|----------------|---|
| 0.005 / 0.025 | βAla, B1, B6-OH, B6-Oxo, NAM, GABA, Sarcosine |
| 0.005 / 0.050 | Cystine |
| 0.005 / 0.075 | Gln, Glu |
| 0.005 / 0.100 | All other analytes |

We're here to help!



+1.844.927.3035
Call us 24/7/365




TechSupport@repligen.com



www.repligen.com/techsupport

For further questions and best practices, please contact your
Repligen field applications scientist.



Thank you

Technical Support | TechSupport@repligen.com

www.repligen.com

