

System Start Up From Long Term Shutdown

PATsmart™ REBEL® System

Support Document

Follow these simple steps when bringing the PATsmart™ REBEL® System back online after a long term shutdown. If any issues are encountered during the steps below, please reach out to our service team at TechSupport@repligen.com or 1-888-927-3035 for assistance.

1. Replace your REBEL consumables with all new - confirm none of the components are expired.
2. Perform a Self-Test in the Settings -> Maintenance Section.
3. Perform a Quantitative Calibration.
4. Set up your batch and run.

Customer Service

Repligen Corporation
685 Route 202/206
Bridgewater, NJ 08807

analytics-support@repligen.com

(908) 707-1009

repligen.com

© 2025 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners.

