

## Preparing your Laboratory for PATsmart™ ZipChip® and Autosampler Installation

This checklist should be completed by the customer to confirm that the on-site requirements for a successful installation of the PATsmart™ ZipChip® Interface and optional Autosampler, have been satisfied. Please return the completed form to your Repligen contact to avoid delays coordinating the installation of your ZipChip Interface. Feel free to email any questions related to this form to [TechSupport@repligen.com](mailto:TechSupport@repligen.com)

### IDENTIFICATION

Name of company / site

Address of site

Name of primary contact

Phone number of primary contact

Date of installation

Primary Application

Secondary Application


### MASS SPECTROMETER MODEL

	Available in Lab	Available for Day of Installation
Thermo Exactive Plus	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive Plus	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive HF	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive HF-X	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive Biopharma	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive HF Biopharma	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive HF-X Biopharma	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Exactive Plus EMR	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive UHMR	<input type="checkbox"/>	<input type="checkbox"/>
Thermo LTQ Orbitrap XL	<input type="checkbox"/>	<input type="checkbox"/>
Thermo LTQ XL	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Exactive Plus	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive	<input type="checkbox"/>	<input type="checkbox"/>
Thermo Q Exactive Plus	<input type="checkbox"/>	<input type="checkbox"/>

## CALIBRATION OF MASS SPEC

Routine calibration of mass spectrometers is required for optimal performance. Please confirm that the calibration of the system has been performed within 48 hours of the installation date. One of the following calibration solutions from Thermo Scientific should be run depending on the mass spectrometer model to be used:

Thermo Scientific Pierce™ LTQ ESI Positive Ion Calibration Solution – Prod #88322 Thermo Scientific Pierce™ LTQ Velos ESI Positive Ion Calibration Solution – Prod #88323 Please refer to the ThermoFisher Scientific website for more information on the calibration solutions. The following spectra should have been generated:



## INSTRUMENT DIMENSIONS AND SPECIFICATIONS

### ZipChip Interface

<b>Mounting Configuration</b>	All Thermo Fisher Scientific Exactive, Q Exactive Orbitrap, and LTQ-Orbitrap MS instruments.
<b>Software</b>	ZipChip software for configuration of injection parameters. Integration with Thermo Fisher Scientific Xcalibur software for data collection, processing and reporting.
<b>Dimensions</b>	8 in (h) x 8 in (w) x 11 in (d)
<b>Weight</b>	11.8 lbs
<b>Power Requirements</b>	110/240V/75W

### Autosampler

<b>Software</b>	ZipChip software for configuration and set up of sequences and run parameters. Routines for automated chip priming and unattended operation.
<b>Dimensions</b>	14 in (h) x 12 in (w) x 22 in (d)
<b>Weight</b>	46.3 lbs
<b>Stackable weight</b>	143 lbs
<b>Power Requirements</b>	110/240V/200W

## CONSUMABLES AND REAGENTS

### ZipChip Qualification Kit - Part Number 850-00043

The kit should be ordered with the instrument, if this does not occur please order 2 weeks prior to install. The kit contains the following:

- 1 set of metabolite reagents
- 2 HS Chip
- Checkout Standards

## CONTACTS

Please contact us if you need any help or support prior to the install:

TECH SUPPORT: [TechSupport@repligen.com](mailto:TechSupport@repligen.com) | 1-888-927-3035

SALES SUPPORT: [analytics-sales@repligen.com](mailto:analytics-sales@repligen.com) | 1-857-254-1500

**Customer Signature:** \_\_\_\_\_