

# LONG<sup>R3</sup> IGF-I Boosts IgG Production in CHO

Shashi Kudugunti, W. Roy Lin and James Rusche  
 Repligen Corporation, R&D, Waltham, MA 02453, USA

## Summary

LONG<sup>R3</sup> IGF-I is a human IGF-1 analog containing a 13 amino acid N-terminal extension and substitution of an arginine for glutamic acid at position 3, which results in >1,000 fold reduced affinity for IGF-BPs, thus enhancing its bioavailability and effectiveness in comparison to native IGF-I. LONG<sup>R3</sup>IGF-I binds to and activates Type I IGF-1 receptor, leading to enhanced cell growth and productivity and it is also a common growth factor supplement used in CHO media as a replacement to insulin. LONG<sup>R3</sup> IGF-I is being used in several commercial antibody manufacturing processes.

An industrially relevant mammalian CHO DP12 cell line (ATCC# CRL-12445<sup>TM</sup>) and cell culture conditions (media + feeds) were selected to evaluate and generate data on the efficacy of LONG<sup>R3</sup> IGF-I. The effect of LONG<sup>R3</sup> IGF-I on the productivity of antibody was compared between LONG<sup>R3</sup> IGF-I adapted and Non- LONG<sup>R3</sup> IGF-I adapted CHO DP12 cells. The fed-batch studies were carried out in shake flasks and results suggested that LONG<sup>R3</sup> IGF-I adapted cells produced two fold increase in the protein production than Non- LONG<sup>R3</sup> IGF-I adapted cells. In addition, the benefit to cost ratio was demonstrated to be high.

LONG<sup>R3</sup> IGF-I can be used to enhance productivity while maintaining a serum free chemically defined media formulation. This provides a robust and cost effective production process.

## Method

### Fed-Batch Method

A frozen vial from LONG<sup>R3</sup> IGF-I adapted CHO DP12 cells was thawed into shake flask containing CD OptiCHO medium with 100 ng/mL LONG<sup>R3</sup> IGF-I, 200nM Methotrexate and 4mM Glutamax. After 4-7 days (at cell density 5-8 e6cells/mL), the cells were inoculated into fed-batch shake flasks containing CD OptiCHO medium with 100 ng/mL LONG<sup>R3</sup> IGF-I, and 4mM Glutamax. A similar thawing method was used for the Non- LONG<sup>R3</sup> IGF-I adapted CHO DP12 cell line vial.

Feeding strategies were similar for both LONG<sup>R3</sup> IGF-I adapted and Non- LONG<sup>R3</sup> IGF-I adapted CHO DP12 cells as listed in the below table.

Fed-Batch Conditions	
Cell Lines	LONG <sup>R3</sup> IGF-I Adapted DP12 Cells Non-LONG <sup>R3</sup> IGF-I Adapted DP12 Cells
Media	CD OptiCHO <sup>TM</sup> Medium (Thermo Fisher# 12681029)
Media supplements	With and Without LONG <sup>R3</sup> IGF-I (100 ng/mL) Methotrexate (200 nM)* Glutamax (4 mM)
Feeds	Day3: Feed A (5%), + Feed B (5%) Day6: Feed A (5%) + Feed B (5%) Day9: Feed A (5%) + Feed B (5%) Day12: Feed A (5%) + Feed B (5%) Feed A (Thermo Fisher# A10234-01) Feed B (Thermo Fisher# A10240-01)
Glucose Feeds	No additional glucose feeds (Feeds already contain high concentrations)
Seeding Density	0.4 E6 Cells/mL
Antibody	Human Anti-IL-8

Note\*: Methotrexate was used only for seed expansion and not in the fed-batch shake flasks

### Fed-Batch Conditions

Flask No.	Fed-Batch Shake Flasks Description	LONG <sup>R3</sup> IGF-I
1	Non-LONG <sup>R3</sup> IGF-I adapted cells (-LONG <sup>R3</sup> IGF-I)	No
2	Non-LONG <sup>R3</sup> IGF-I adapted cells (+LONG <sup>R3</sup> IGF-I)*	Yes
3	LONG <sup>R3</sup> IGF-I adapted cells (+LONG <sup>R3</sup> IGF-I)	Yes


Note\*: In shake flask 2, LONG<sup>R3</sup> IGF-I was added to fed-batch culture on Day 0 in spite of the cells being Non-LONG<sup>R3</sup> IGF-I adapted

## Conclusions


- LONG<sup>R3</sup> IGF-I adapted CHO DP12 cells produced 2-fold increase in the protein yield than Non- LONG<sup>R3</sup> IGF-I adapted cells.
- LONG<sup>R3</sup> IGF-I adapted CHO cells produced 100% increase in IgG production with less than 10% in costs, when compared to the Non- LONG<sup>R3</sup> IGF-I adapted cells.
- Increased production of IgG was observed only after adapting the cells to LONG<sup>R3</sup> IGF-I. Supplementing the LONG<sup>R3</sup> IGF-I to the media without cell line adaption might not be helpful.
- LONG<sup>R3</sup> IGF-I adapted CHO DP12 cells did not lose cell specific production rate in a 25 day ATF perfusion culture.

## Results

### Viable Cell Density (VCD) and Viability





### Protein Production




### Cell Specific Production Rate of LONG<sup>R3</sup> IGF-I Adapted Cells

This perfusion experiment was conducted separately with LONG<sup>R3</sup> IGF-I adapted DP12 cells inoculated into a bioreactor connected to Repligen ATF2 perfusion system to determine cell specific production rate. The run was conducted without methotrexate and VCD was maintained at 40-60 E6 cells/mL.



## Cost Benefit Analysis



Note: Graph was prepared accounting for approximate fed-batch expenses on a large scale production (Media: 10\$/L and LR3 20\$/mg)