

# MabSelect SuRe™ ELISA Kit Performance Summary

## 1. Introduction

The Protein A ELISA Kit (catalog number 9333-1) is intended for the detection and quantitation of residual MabSelect SuRe™ ligand. Our Protein A detection ELISA Kit has been developed for those customers who require a highly sensitive assay to measure small amounts of contaminating MabSelect SuRe™ ligand in antibody products. Testing for contamination occurs in several different phases of development and commercial manufacturing that may include:

- Process development for leaching characteristics of the resin under specific conditions
- Manufacturing, typically from eluted samples taken throughout several points in the purification process
- Finish product release to document process containment levels and lot-to-lot consistency

The following summary report contains performance data collected from the evaluation of the Protein A ELISA Kit in the presence of human Immunoglobulin G (hIgG). The data presented here demonstrates the Protein A ELISA Kit's:

- Ability to detect the MabSelect SuRe™ ligand in the presence of up to 0.5 mg/ml hIgG in a PBS-T buffer
- Percent recovery (accuracy), inter and intra assay precision, limit of quantitation and limit of detection

## 2. Results Summary

### MabSelect SuRe™ in the Presence of hIgG

The performance of the Protein A ELISA Kit was evaluated when detecting the MabSelect SuRe™ ligand in the presence of hIgG compared to a standard containing no hIgG. All spiked samples had a final hIgG concentration of 0.125 mg/ml (following final dilution into the assay plate). Each sample was prepared in replicates of 8 and three separate ELISAs were performed according to the kit's standard protocol.

### Data Handling

Standard curve data points were fitted to a 4-parameter fit analysis. This equation allowed back-calculation of sample rPA concentrations and calculation of LoQ values. Percent recovery was calculated as follows:

$$\% \text{ Recovery} = \frac{\text{Calculated Conc}}{\text{Theoretical Conc}} \times 100$$

### Intra-Assay Precision

Table 1. Intra-Assay Precision for Standard Curve samples

Conc (ng/ml)	Avg % CV
1.6	5.1
0.8	2.3
0.4	3.3
0.2	4.0
0.1	5.2
0.05	7.5

MabSelect SuRe™ is a trademark of GE Healthcare.

# MabSelect SuRe™ ELISA Kit Performance Summary

Table 2. Intra-Assay Precision for samples containing hIgG

Conc (ng/ml)	Calculated Conc	Avg %CV
1.2	1.25	6.1
1.0	1.01	5.3
0.8	0.76	4.1
0.6	0.64	4.9
0.4	0.37	5.4
0.2	0.18	7.7
0.1	0.09	11.0
0.05	0.05	18.0

## Inter-Assay Precision

Table 3. Inter-Assay Precision for Standard Curve samples

Conc (ng/ml)	Avg % CV
1.6	0.1
0.8	0.2
0.4	0.8
0.2	1.6
0.1	2.1
0.05	2.7

Table 4. Inter-Assay Precision for samples containing hIgG

Conc (ng/ml)	Calculated Conc	Avg %CV
1.2	1.25	7.6
1.0	1.01	0.9
0.8	0.76	2.6
0.6	0.64	8.6
0.4	0.37	6.3
0.2	0.18	13.0
0.1	0.09	10.4
0.05	0.05	9.0

## Accuracy

Table 5. Accuracy for Standard Curve samples

Conc (ng/ml)	Avg % Error	Avg % Recovery
1.6	0.0	100
0.8	0.4	100
0.4	-0.7	99
0.2	-1.3	99
0.1	2.1	102
0.05	7.3	107

MabSelect SuRe™ is a trademark of GE Healthcare.

# MabSelect SuRe™ ELISA Kit Performance Summary

Table 6. Accuracy for samples containing hIgG


Conc (ng/ml)	Calculated Conc	Avg % Error	Avg % Recovery
1.2	1.25	3.8	104
1.0	1.01	1.2	101
0.8	0.76	-4.5	96
0.6	0.64	5.9	106
0.4	0.37	-8.2	92
0.2	0.18	-8.0	92
0.1	0.09	-8.8	91
0.05	0.05	-7.8	92

## Limit of Quantitation (LoQ)

The LoQ determined from the standard curve was calculated to be 0.025 ng/ml. The LOQ for the MabSelect SuRe™ ligand spiked samples in the presence of hIgG was 0.115 ng/ml or 0.920 ng/mg (0.92 ppm). Based on this data set the sensitivity of the kit when detecting MabSelect SuRe™ ligand in the presence of immunoglobulin is < 1 part per million (ppm).

## Linearity of Dilution

Figure 1. Linearity of dilution graph showing high correlation of accuracy throughout the range of MabSelect SuRe™ ligand concentrations tested.



MabSelect SuRe™ is a trademark of GE Healthcare.

# MabSelect SuRe™ ELISA Kit Performance Summary

## Conclusions

Intra-assay data indicated that precision of the three assays was within acceptable parameters, with % CV values of the standards ranging from 2.3-7.5 %. This set of experiments had a standard LoQ of 0.025 ng/ml.

Intra-assay data indicated that samples containing hIgG were assayed with relative accuracy and precision, with % CV values from 11.0-4.1%. The LoQ calculated for the samples was 0.115 ng/ml. Accuracy, gauged by % relative error, was from 7.3 to -1.3 for standards, and from -8.8 to 5.9 for samples.

Inter-assay data had similar findings. The % CV between standard curve samples was less than 2.7 for all concentrations, indication of precision across the three assays. The samples containing hIgG had a high degree of precision as well with % CV values in the range of 13.0-0.9% for concentrations above the LoQ.

## 3. Explanation of Calculations

### Precision (%CV)

Precision was calculated by determining the standard deviation between rPA spiked sample data points and dividing by the mean value. According to the 'Guidance for Industry: Bio-analytical Method Validation' text, precision should be within 15%.

### Intra-Assay Precision

The intra-assay precision was calculated for each rPA spiked sample concentration by averaging the %CV values across all assays.

### Inter-Assay Precision

The inter-assay precision was calculated for each concentration point by determining the standard deviation between calculated results from each of the three assays, then dividing by the mean value.

### Limit of Quantitation (LoQ)

The limit of quantitation (LoQ) was defined as 10 times the standard deviation of 0 ng/ml sample. The standard deviation of the 0 ng/ml OD value was multiplied by 10 then added to base 0 ng/ml OD value. The LoQ was then generated by entering the summed value into the standard curve 4-parameter fit equation. For each kit the LoQ was reported as ng Protein A per ml (ng/ml) buffer, and ng Protein A per mg hIgG (ppm) for rPA spiked samples run in presence of hIgG.

### Limit of Detection (LoD)

The limit of detection (LoD) was defined as 3 times the standard deviation of 0 ng/ml Protein A sample. The standard deviation of the 0 ng/ml OD value was multiplied by 3 then added to base 0 ng/ml OD value. The LoD was then generated by entering the summed value into the standard curve 4-parameter fit equation. For each kit a LoD was reported as ng Protein A per ml (ng/ml) buffer and in parts per million (ppm).

### Accuracy

Accuracy is described as the % recovery determined by the assay compared to the theoretical spiked concentration.

MabSelect SuRe™ is a trademark of GE Healthcare.