

Automated Parallel Chromatographic Separations in Process Development

Tim Schroeder, Juergen Friedle

Atoll GmbH, Ettishofer Straße 10, D-88250 Weingarten, Germany

Summary

Development and optimization of chromatographic separation conditions is normally a time consuming task in the large scale production of proteins such as monoclonal antibodies. In this study a 96 well formatted MediaScout® RoboColumn® array was adapted for automated operation in a modified commercial liquid handling workstation and used for the development of a purification strategy for a recombinant, secretory enzyme by screening different cation exchange chromatography media according to optimize yield and purity (Fig. 1). The quality of the purified enzyme was checked by using cutting-edge high throughput LabChip-SDS Electrophoreses (Caliper).

The combined approach of Atoll's MediaScout® chromatography tools and Tecan's liquid handling workstation Freedom EVO® enables to speed up the development of a purification process by nearly one order of magnitude.





FIGURE 1 API development. Process flow.

Experimental Setup



Patent Pending PCT/DE 2006/000708

FIGURE 2 MediaScout® RoboColumn®.



FIGURE 3 Liquid handling workstation Tecan Freedom EVO. Modified for use with RoboColumn®.





FIGURE 4 a LabChip GXII System. b Microfluidic chip.

Purification

FIGURE 5 Automated small scale purification of enzyme using RoboColumn®. Elution profile.

- Fractogel® EMD SO₃⁻ shows a higher overall protein binding capacity than Capto™ S.
- Elution profile does not include data about the purity of the desired enzyme. → SDS Electrophoreses




Column	RoboColumn®
V _{Column}	200 µl
Media 1	Fractogel® EMD SO ₃ ⁻
Media 2	Capto™ S
V _{Sample}	6.0 ml
Flow Rate	150 cm/h
Binding buffer	10 mM MOPS, pH 6.8 50 mM NaCl
Elution buffer 1	10 mM MOPS, pH 6.8 150 mM NaCl
Elution buffer 2	10 mM MOPS, pH 6.8 400 mM NaCl
Elution buffer 3	10 mM MOPS, pH 6.8 800 mM NaCl

SDS Electrophoreses

FIGURE 6 SDS Electrophoreses after purification of enzyme using LabChip technology.

- Capto™ S shows a more specific binding behavior and therefore higher purity for the desired enzyme.
- Fractogel® EMD SO₃⁻ is binding a lot of undesired proteins (impurities) as well.

Media	Fractogel® SO ₃	Capto™ S
Yield	++	+
Purity	--	++
Σ	0	+++



Conclusions

- Atoll's 96 MediaScout® RoboColumn® array was successfully adapted on Tecan's liquid handling workstation Freedom Evo® for chromatographic applications.
- Small scale automated high throughput separations using bio-chromatography were successfully applied for screening of cation exchange chromatography media in process development.
- The desired enzyme was successfully concentrated and purified using MediaScout® RoboColumn® packed with Capto™ S.
- Using the combined approach of Atoll's MediaScout® RoboColumn®, Tecan's liquid handling workstation Freedom EVO and Caliper's LabChip technology enables to reduce costs per experiment significantly due to saving process time, API, process-relevant products and solvents.