

Intensified Seed Expansion and High Productivity Harvest of Fed-batch Cultures Using the XCell™ ATF System

Shashi Kudugunti*, Jyoti Amatya, Dan Diggins, Amit Dutta, & James Peyser; Repligen Corporation, Waltham, MA 02453, USA
* Presenter


Introduction

A conventional fed-batch (FB) process involves a multi-stage seed expansion and a large footprint of clarification tools to isolate the product of interest from cell culture fluid containing cells and cell debris. The current standard clarification technology of centrifugation and depth filtration often requires open processing and a large footprint. This multi-stage process is cumbersome, difficult to scale, and regularly impacts overall yield. With an increase in demand for production of biotherapeutics, it is essential to simplify the traditional FB process to result in higher yield, shorter process time and smaller footprint.


Repligen developed an intensified FB process using XCell ATF technology. This process involves inoculating the FB bioreactor at high seeding density followed by a novel clarification method, High Productivity Harvest (HPH). This case study presents the N-1 bioreactor in a perfusion mode allowing the FB bioreactor to be inoculated at 20X higher seed density than a traditional process. Increased seed density reduces the FB culture duration by 50%. In addition, implementing the HPH clarification process using the XCell ATF system results in higher VCD and viability due to the replacement of nutrients and removal of toxins. The addition of intensified seed expansion provides a 190% boost in the overall yield in half the process time.

Maintaining high viability throughout the process provides a healthier environment for improved product quality. The entire process is performed as a single step in an integrated, enclosed, bioburden free environment, yielding clarified product ready for downstream purification.

Conventional vs. Intensified Fed-batch Process



Novel High Productivity Harvest Process using XCell™ ATF



Materials & Methods


High expressing mammalian CHO-GS, was selected to evaluate the HPH process. These cells were grown in Ex-CELL FB medium supplemented with 100 ng/mL LONG®R³ IGF-I, 4g/L Poloxamer, and 4 mM Glutamax. This cell line is reported to express IgG, trastuzumab. All cultures were conducted using 3L (w.v 1.2L) Applikon glass bioreactors. The cell line and the method for traditional FB process (including feeding strategy) is provided by ATUM¹ & Horizon Discovery².

Abbreviations: DD Clarification (Drawdown Clarification), FB (Fed-Batch), IGF-1 (Insulin like Growth Factor 1), HPH (High Productivity Harvest), HCP (Host Cell Protein), LMH (L/m²*H), SD (Seeding Density), VCD (Viable Cell Density), VVD (Vessel Volume/day)


Components /Conditions	FB Traditional Low SD	FB_HPH Low SD	FB Traditional High SD	FB_HPH High SD
(N-1) Process	(N-1) Batch Duration: 4 days Achieved VCD: 6E6 cells/mL	(N-1) Perfusion Duration: 11 days Achieved VCD: 100E6 cells/mL	(N-1) Batch Duration: 4 days Achieved VCD: 6E6 cells/mL	(N-1) Perfusion Duration: 11 days Achieved VCD: 100E6 cells/mL
FB Seeding Density	0.5E6 cells/mL	10E6 cells/mL	0.5E6 cells/mL	10E6 cells/mL
FB Duration	14 days	7 days	14 days	7 days
XCell™ ATF System	N/A	XCell™ ATF2 PES 0.2µm	N/A	XCell™ ATF2 PES 0.2µm
Feed/Perfusion Rate	7.5% Feed: Day 3, 5, 7, 9, 11 & 13 ↓ Day 14: Clarification	7.5% Feed: Day 3, 5, 7, 9, 11 & 13 ↓ Day 10: 1 VVD HPH Day 14: DD Clarification	7.5% Feed: Day 0, 1, 2, 3, 4, 5, & 6 ↓ Day 7: Clarification	7.5% Feed: Day 0, 1, 2, 3, 4, 5, & 6 ↓ Day 3: 1 VVD HPH Day 7: DD Clarification
Drawdown (Flux)	N/A	10LMH	N/A	10LMH
Shear Rate	N/A	2000s ⁻¹	N/A	2000s ⁻¹

Results


Viable Cell Density & Viability




Cell Specific Protein Productivity



Total Accumulated Protein



Product Quality



Conclusions

A production bioreactor seeded with a high cell density and operated with the HPH process results in significant increase in product yield, improved process efficiency and smaller operational footprint. The productivity and process benefits include:

- Up to 190% boost in protein production
- 50% reduction in production culture duration
- Eliminate centrifugation and depth filtration equipment and operations
- Healthier culture translates to lower HCP and lower process related impurities
- Closed system and single-step process
- Harvested material is 0.2µm filtered and ready for purification