

Intentionally designed for high-pressure applications

The tubing and hose clamp meets the requirements for critical fluid applications within biotech, pharmaceutical, food, and other high-purity applications. User-friendly with quick swing-away lock-down, it can be installed on existing assemblies without disturbing the process fluid path.



- **Lightweight and reusable**
- **Autoclavable** to 135° C
- **Gamma stable** to 50 kGy
- **UPS Class VI material**
- **Swing-away locknut**
- **Wide range** of tubing sizes
- **Validated designs**

