

Reusable TFF cassettes for ultrafiltration/diafiltration

Delivering equivalent performance, scalable and interchangeable with TangenX[®] SIUS[®] Single-use Cassettes

TangenX[®] PRO Reusable Tangential Flow Filtration (TFF) Cassettes may be cleaned and reused multiple times while maintaining reproducible cassette performance. Validated for GMP use, TangenX[®] PRO Cassettes come in two industry standard format footprints with accompanying holders. They are available using ProStream or HyStream membranes, with choice of pore size and flow-channel configuration.



For downstream concentration and diafiltration

Interchangeable with TangenX[®] SIUS[®] Single-use Cassettes

Up to 30% less cost than other reusable cassettes

Scalable from process development to commercial scale

Compatible with all industry holders



TangenX[®] PRO PD Cassettes for process development

Membrane surface areas
0.01 m² 0.02 m² 0.10 m²

TangenX[®] PRO Cassettes for production scale

Membrane surface areas
0.5 m² 1.5 m² 2.5 m²

Revised 6/11/2019

