

REBEL SMARt® Kit

# Cell Culture Media Analysis at Your Fingertips

With REBEL SMARt® Kit\* you can quickly analyze key components in cell culture media in about 15 minutes per sample, providing actionable insights for accelerated process development.

**NEW!**

Extended shelf life,  
focused analyte  
panel

## All-Inclusive Kit:

- SMARt Standards (3 strips)
- REBEL® XT Chip (1)
- REBEL® Diluent (2)
- REBEL BGE (2)
- 96-Well Plate & Covers (2)



**Contact us**

\* For use with the REBEL® XT Media Analyzer only

## Focused Analyte Panel: 19 amino acids, 2 dimers, and choline

- |                    |                 |                 |              |
|--------------------|-----------------|-----------------|--------------|
| • Alanine          | • Cystine       | • Leucine       | • Threonine  |
| • Alanyl-Glutamine | • Glutamic acid | • Lysine        | • Tryptophan |
| • Arginine         | • Glutamine     | • Methionine    | • Tyrosine   |
| • Asparagine       | • Glycine       | • Phenylalanine | • Valine     |
| • Aspartic acid    | • Histidine     | • Proline       |              |
| • Choline          | • Isoleucine    | • Serine        |              |

## Flexible Ordering. Reliable Delivery. Just What Your Lab Needs.

Choose the purchasing and shipping option that works best for you:

### One-Time Purchase:

Order kits as needed at full price—  
quick and easy.

### Subscription Savings:

Order 6, 12, or more kits per shipment at a  
discounted rate, delivered on your schedule.

**repligen.com**

© 2025 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners. | DOC0429 eRev. 1.0 21 May 2025

