

PATsmart™ REBEL® & REBEL® XT Consumables Kits

Storage and Expiration Specifications

Specification Sheet

This guide serves as a reference for proper storage conditions and shelf life for chemical components received in consumables kits for the PATsmart™ REBEL® and REBEL® XT Systems. **Please note:** The REBEL standards are shipped separate from the remainder of the kit in order to maintain cold storage (2°C–8°C) in transit. The REBEL standards pouch should be refrigerated immediately upon arrival. The contents of the REBEL kit box do not require refrigeration.

Part Number	Description	Components	Shelf Life*	Storage Temperature	Working Life**
850-00050	REBEL SMA Kit	1 x REBEL chip 2 x 250 mL BGE 2 x 125 mL Diluent with internal standard 2 x 96-well microplates SMA Standards	BGE & Diluent: 13 months Standards: 4 months	BGE and diluent: room temperature Standards: 2–8 °C. Do not freeze.	BGE and diluent: 30 days Standards: 14 days
850-00083	REBEL SMA Kit without Standards	1 x REBEL chip 2 x 250 mL BGE 2 x 125 mL Diluent with internal standard 2 x 96-well microplates	BGE & Diluent: 13 months	BGE and diluent: room temperature	BGE and diluent: 30 days
850-00059	REBEL IQOQ Kit	1 x REBEL chip 1 x 250 mL BGE 1 x IQ/OQ Standard Kit (2 vials Performance Qualification std; 5 vials Quantitation std; 5 vials IQ/OQ sample)	BGE: 13 months Standards: 4 months	BGE: room temperature Standards: 2–8 °C. Do not freeze.	BGE: 30 days Standards: 14 days
850-00064	REBEL SMA Standards	3 x 5 vials Quantitation standard 3 x 2 vials Performance Qualification Standard	4 months	Shipped cold and must be refrigerated (2–8 °C) upon arrival. Do not freeze.	14 days
850-00092	REBEL SMARt Kit	1 x REBEL XT chip 2 x 250 mL BGE 2 x 125 mL Diluent with internal standard 2 x 96-well microplates SMARt Standards	BGE & Diluent: 13 months Standards: 9 months	BGE and diluent: room temperature Standards: 2–8 °C. Do not freeze.	BGE and diluent: 30 days Standards: 14 days
850-00118	REBEL XT IQOQ Kit	1 x REBEL XT chip 1 x 250 mL BGE 1 x IQ/OQ Standard Kit (2 vials Performance Qualification std; 5 vials Quantitation std; 5 vials IQ/OQ sample)	BGE: 13 months Standards: 9 months	BGE: room temperature Standards: 2–8°C. Do not freeze.	BGE: 30 days Standards: 14 days
850-00094	REBEL SMARt Standards	3 x 5 vials Quantitation standard 3 x 2 vials Performance Qualification Standard	9 months	Shipped cold and must be refrigerated (2–8 °C) upon arrival. Do not freeze.	14 days

*Shelf life from manufacture date.

**Working Life: Maximum recommended use from initial opening when capped/sealed and stored properly

repligen.com

REBEL AND REBEL XT ARE FOR RESEARCH USE ONLY

REBEL and REBEL XT are subject to export controls including those of the Export Administration Regulations of the U.S. Department of Commerce, which may restrict or require licenses for the export of product from the United States and their re-export to and from other countries. © 2026 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners. | DOC0585 eRev. 1.0 19 Aug 2026

