

# Vial and Plate Recommendations

## PATsmart™ REBEL® System

Support Document

Below are the approved consumables that have been tested and certified to work with the PATsmart™ REBEL® System.

### 96 Well Plates

- [249944](#) - Thermo Scientific, Nunc 96-Well Polypropylene Sample Processing & Storage Microplates (non-sterile)
- [249946](#) - Thermo Scientific, Nunc 96-Well Polypropylene Sample Processing & Storage Microplates (sterile)

### Well Plate Covers

- [AB0626](#) - Thermo Scientific, Adhesive PCR Plate Foils

### Autosampler Vials

#### Amber Autosampler Vial Kits

- [6AK594W](#) - Thermo Scientific, SureSTART Specification Certified Screw Vial and Cap Kits, Level 2 High-throughput Applications, Pack of 100, Amber Glass

#### Amber Individual Vial Components

- [6ASV9-2P](#) - Thermo Scientific, SureSTART 2 mL Glass Screw Top Vials, Level 2 High-Throughput Applications (Amber Glass - Vials only)
- [6ASC9ST1](#) - Thermo Scientific, SureSTART 9 mm Screw Caps, Level 2 High-throughput Applications (Blue Caps Only)
- [6ASC9ST1R](#) - Thermo Scientific, SureSTART 9 mm Screw Caps, Level 2 High-throughput Applications (Red Caps Only)

*\*Clear/Colorless autosampler vials do not provide protection for photosensitive analytes. If samples are not prepared and stored in the proper environment, photodegradation of some analytes may occur. Analytes that may experience photosensitivity include, but may not be limited to: Cystine, His, Lys, Met, Phe, Thiamine (B1), Trp, and Tyr.*

### Customer Service

Repligen Corporation  
685 Route 202/206  
Bridgewater, NJ 08807

[TechSupport@repligen.com](mailto:TechSupport@repligen.com)

(908) 707-1009

[repligen.com](https://repligen.com)

© 2025 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners.

