

PATsmart™ ZipChip® Software PC Requirements

Guidelines

The following are the recommended minimum PC requirements for running the PATsmart™ ZipChip® software:

- Windows 7 or later
- 64-bit Operating system
- Dual Core CPU, at least 2 GHz
- 1 GB free disk space on hard drive
- 8 GB RAM

Customer Service

Repligen Corporation
685 Route 202/206
Bridgewater, NJ 08807

analytics-support@repligen.com

(908) 707-1009

repligen.com

© 2025 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners.

