


Specifications



Model SCL-0020-SCLR	
Maximum capacity	20 kg
Readability	0.1 g
Tare range	-20 kg
Repeatability	0.1 g
Linearity	±0.3 g
Stabilization time	1 s
Adjustment	External
Display	LCD (with backlight)
Keypad	5 keys
Protection class	IP 43
USB-A	1
USB-B	1
RS 232	2
Power consumption	1.32 VA
Power supply	Pump (direct connection with octocable)
Weighing pan dimensions	7.7 in x 7.7 in
Weighing device dimensions	13 in x 8 in x 3.8 in
Net weight	6 lb
Gross weight	9.5 lb
Packaging dimensions	15 in x 21 in x 17 in

Dimensions



Revised 04-14-2014

