


Easily make sanitary connections

The Lightweight Sanitary Clamp, ideally designed for Single Use systems in the Bioprocess, Pharmaceutical, and Cell and Gene Therapy Industries



- **High-purity** TC Clamp
- **User-friendly**
- **Hygienic**
- **Compatible** with ASME-BPE flange and gasket standards
- **Suitable** for non-metallic connections
- **Available sizes:** 3/4 in , 1-1/2 in, 3 in,



Secure bottles and protect final product

Drug substance security

Easily protect drug substance and add a tamper-proof seal to any bottle, container or carboy used in biopharmaceutical production. The simple, easy-to-use lockout cover snaps in place within seconds, adding another level of risk mitigation.

Avoiding issues with bottles

Current storage bottles are convenient, but simple screw top designs can easily and inadvertently be tampered. Common bottle failures occur at the bottle cap tubing system. The Bottle Locker protects tubing closures and eliminates accidental exposure of bottle contents to the environment.



Applications

- Freeze Thaw
- Drug Product Protection
- Viral Containment
- Product Transportation
- Cell and Gene Therapy



Mitigate the risk of shipping bottles

with common screw caps

Next level product protection

The Bottle Locker was designed as a tamper evident system for active ingredients that aligns with similar industry regulations currently required by law on final product.

Features

- Fits Industry standard bottles
- Supports 100mL to 20 Liters

Ordering instructions

| | |
|-----------------------------|--------------------------------------|
| 1. Select bottle type | Media bottle, BioTainer, Carboy |
| 2. Confirm closure thread | 38 - 430, 48 mm, GL 45 |
| 3. Determine closure height | Standard cap, cap with tubing insert |
| 4. Add custom requirements | Material of construction, lockout |

